MerFill Bond SE

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info@famedical.net

www.famedical.net

MerFill Bond SE is a self-etching single-component light curing adhesive which provides excellent bonding to dentine, enamel, composites, metal and amalgam.

MerFill Bond SE contains phosphoric acid methacrylate esters and maleic acid-modified methacrylates which being placed in aqueous surroundings will provide the moderate etching the surfaces of dentin and enamel to allow for the chemical and micromechanical bonding of a restorative material to tooth tissue. These chemically active components can react with tooth hydroxyapatite by forming a chemical bond with the calcium ions. In process of bonding of MerFill Bond SE to dentin the dissolving of inorganic smear layer and the demineralization of the intertubular dentin tissue are taken place. In the same time the "polymerizable" adhesive component penetrates into tubular structure of demineralized dentin for forming the resin tags which provide the strong adhesion between tooth tissue and restoration.

•• MAIN ADVANTAGES

- Guarantees a simple, fast and economical application technique.
- Contains components that provide elasticity when working with all types of composites by creating an elastic hybrid layer. Such elastic hybrid layer reduces the risk of microleakage due to high elasticity and resistance to kinking. It also contains surfactants (surface active co-polymerized components) that facilitate penetration of the resin into the dentinal tubules.
- A special co-polymerized solvent in combination with a very small amount of ethyl alcohol as a carrier solves the problem of rapid evaporation of the excess (undesirable) solvent and at the same time ensures the penetration of the adhesive into the dentinal tubules.
- Does not contain acetone or other compounds that irritate the pulp.

•• INDICATIONS

- Self-etching adhesive system for restoring all classes of cavities (according to Black) with the use of composite and compomer restorative materials (including flowing materials).
- Adhesive system for fixing indirect restorations with composite and compomer materials and ceramics with the use of composite double-curing cements.
- Protective varnish for the cervical areas of the teeth.

•• CONTRAINDICATIONS

- Using for fixing restorations with composite or compomer materials of double or chemical curing.
- Direct or indirect coating of pulp without the use of lining materials.
- Use for treatment of patients with an allergic reaction to methacrylates.

•• MODE OF APPLICATION

• INSULATION

It is recommended to use a cofferdam for insulation.

• TOOTH PREPARATION

Make the tooth preparation using standard methods. Carefully remove damaged tooth tissue. Use water to remove debris formed after preparation, dry with air or cotton pellets. Do not overdry. The surface should be slightly glossy.

• ETCHING

To increase adhesion when working on the enamel area of the tooth, additional etching of enamel areas is recommended before applying the adhesive.

Unprepared enamel must be thoroughly cleaned with a polishing paste or with the intraoral sandblasting machine. Etching of the enamel is fulfilled with the help of standard etching gels for 15-20 seconds.

The prepared enamel can be selectively etched to increase adhesion. Apply etching gel to the enamel, including the overhanging edges of the enamel. Leave for 15-20 seconds. Rinse well with water, then air dry. Do not overdry.

• PULP PROTECTION

To avoid irritation of the pulp, cover the area near the pulp with calcium-containing pads.

• ADHESIVE APPLICATION AND RESTORATION PROCEDURE.

Apply MerFill Bond SE with a brush rubbing the material well into the prepared cavity.

Wait 20-25 seconds. Dry with air for 5 seconds. Do not rinse. Light-cure for 10-20 seconds. Then, perform the restoration with light-curing composite materials- either directly or indirectly. Follow the instructions of the material manufacturer.

When fixing indirect restorations to composite cement, prepare the restoration for fixation in accordance with the requirements of the manufacturer of the material for indirect restorations.

• SPECIAL RECOMMENDATIONS

When restoring V class cavities (cervical defects) the adjacent edges of the enamel should be bevelled and etched in the usual way. Before etching dentin and adjacent enamel must be cleaned with a pumice-water suspension (free of oil traces), washed and lightly dried.

• APPLICATION TO REDUCE THE HYPERSENSITIVITY OF THE CERVICAL AREAS

Clean the cervical area of the tooth with increased sensitivity, using a preventive paste and a rubber cup. Rinse the paste with water. Dry with a jet of air without oil. Avoid drying the dentine. Apply MerFill Bond SE in the same way as in the case of direct restorations.

•• WARNING

- Eugenol-containing materials should not be used with MerFill Bond SE, since eugenol can adversely affect the polymerization reaction.
- Contact with soft tissues may cause allergic contact dermatitis.

•• PRECAUTIONARY MEASURES

1. During the work with the material use special glasses, gloves, clothes and mask. Goggles are recommended for the patients also.
2. Observe the security measures provided during the work with powdery, siliceous substances. Use goggles, a mask and gloves.

•• WARRANTY

The manufacturer warrants the quality of manufactured products. The adverse events inflicted by violation of user manual, storage conditions and other events inflicted by non-stipulated usage of the material are not the subjects of warranty. The customer is responsible for determination of suitability of this product for user's application. Warranty conditions: the product does not comply with requirements declared by manufacturer. In this case the manufacturer replaces the defective material within warranty period.

•• CONTRAINDICATIONS

Allergic reactions in the anamnesis of the patient especially on methacrylic resins or on any other component of dental materials.

•• LIMITATION OF LIABILITY

The manufacturer's liability is limited by only cases stipulated by direct legislation of the country.

•• COLLATERAL REACTIONS

The product can cause irritation of eyes, skin and mucous membranes.

For professional use in dentistry only!

•• STORAGE

- Keep syringe or bottle with MerFill Bond SE tightly closed after each using. Don't allow direct sunlight on the product.
- The shelf-life of MerFill Bond SE is two years under the temperature 2 – 8 °C (32 – 46°F).
- The shelf-life of MerFill Bond SE is about 8 months if store temperature is 23 °C (73 °F).
- DO NOT FREEZE!

•• RECYCLING

Dispose of the medical device in accordance with local / regional / national / international legal requirements.

•• PACKING

- MerFill Bond SE in syringes of 3 ml or in bottles 6 ml.

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